

Creatine Kinase (CK) Activity Assay Kit (Spectrophotometry)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Creatine Kinase (CK) (EC 2.7.3.2), also known as creatine phosphokinase, is primarily found in tissues such as the heart, muscles, and brain. It reversibly catalyzes the transphosphorylation reaction between creatine and ATP, it is an important kinase directly involved in cellular energy metabolism, muscle contraction, and ATP regeneration.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481