

## Citrate synthase Assay Kit (Microanalysis)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years  
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Citrate Synthase (CS) (EC 2.3.3.1) is widely present in the mitochondrial matrix of animals, plants, microorganisms, and cultured cells. It is the first rate-limiting enzyme of the tricarboxylic acid (TCA) cycle and one of the major regulatory sites of the cycle.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481