

Citrate synthase Assay Kit (Spectrophotometry)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Citrate Synthase (CS) (EC 2.3.3.1) is widely present in the mitochondrial matrix of animals, plants, microorganisms, and cultured cells. It is the first rate-limiting enzyme of the tricarboxylic acid (TCA) cycle and one of the major regulatory sites of the cycle.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481