

CSP Assay Kit (Spectrophotometry)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Pectin is a major component of the primary cell wall and middle lamella, primarily consisting of protopectin, pectin methyl ester, and pectic acid. Pectin contains galacturonic acid, lactose, arabinose, glucuronic acid, making it the most abundant polysaccharide in the cell walls of many higher plants. Its unique physical and chemical properties influence the taste and quality of plant-based foods. Pectin molecules are cross-linked via Ca^{2+} bridges, ionic bonds, hydrogen bonds, glycosidic bonds, ester bonds, and phenyl ring couplings. Various extraction methods can yield different forms of pectin, such as water-soluble pectin (WSP), ion-bound pectin (ISP), and covalently bound pectin (CSP).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481