

CYB5 Assay Kit (Visible Spectrophotometry)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Cytochrome P450 enzymes are a group of isoenzymes mainly found in the liver, playing a significant role in the metabolism of xenobiotics, particularly drugs and toxins. Cytochrome P450 and cytochrome b5 (CYB5) are two hemoproteins in the P450 enzyme system, and changes in their ratio are closely related to P450 metabolic activity.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481