

DiamineOxidase Assay Kit (Visible Spectrophotometry)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Diamine oxidase (DAO; EC 1.4.3.6) is widely found in animals (intestinal mucosa, lungs, liver, kidneys, etc.), plants, and microorganisms. It catalyzes the oxidation of polyamines into aldehydes. Its activity is closely related to nucleic acid and protein synthesis and can reflect the integrity and extent of damage to the intestinal mechanical barrier.
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481