

Chitinase Assay Kit (Microanalysis)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Chitin is primarily found in the exoskeletons of crustaceans such as shrimp and crabs, the shells of insects, the organs of mollusks (such as squid cartilage), and the cell walls of fungi. Chitinase (EC 3.2.1.14) catalyzes the hydrolysis of chitin and plays a role in defending against fungal infections, making it a major research focus in the study of antifungal disease resistance.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481