

Fatty Acid Synthase Assay Kit (Microanalysis)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fatty Acid Synthase (FAS) is a key enzyme in fatty acid synthesis, catalyzing the formation of long-chain fatty acids from acetyl-CoA and malonyl-CoA. FAS is widely expressed in various tissues and cells, with particularly high expression in mammalian liver, kidney, brain, lung, mammary gland, and adipose tissue.
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481