

Compound TYD-02596

Chemical Properties

CAS No. : 333425-29-9

Formula:

Molecular Weight:

Store at low temperature

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Compound TYD-02596, with CAS No. 333425-29-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Compound TYD-02596 provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481