

Sulfo Cy3(Et)-N3 DIPEA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3-N3,sulfo Cy3(Et)-N3.DIPEA is a CY3-derived dye. CY3 (Sulfo-Cyanine3) is an amine-reactive cyanine orange dye with an Ex/Em of 550/570 nm, commonly used for fluorescent labelling of proteins and nucleic acids.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481