


## Atropic acid

## Chemical Properties

CAS No. :	492-38-6
Formula:	C <sub>9</sub> H <sub>8</sub> O <sub>2</sub>
Molecular Weight:	148.16
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Atropic acid (Acrylic acid) (Ipratropium EP Impurity D) is an impurity of Ipratropium bromide and Atropine Sulfate.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 25 mg/mL (168.74 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7495 mL	33.7473 mL	67.4946 mL
5 mM	1.3499 mL	6.7495 mL	13.4989 mL
10 mM	0.6749 mL	3.3747 mL	6.7495 mL
50 mM	0.135 mL	0.6749 mL	1.3499 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Lyon RC, et al. Stability profiles of drug products extended beyond labeled expiration dates. J Pharm Sci. 2006 Jul; 95(7):1549-60.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481