


Tiplimotide

Chemical Properties

CAS No. :	178823-49-9
Formula:	C87H142N24O21
Molecular Weight:	1860.21
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Tiplimotide is an altered peptide ligand (APL).
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5376 mL	2.6879 mL	5.3757 mL
5 mM	0.1075 mL	0.5376 mL	1.0751 mL
10 mM	0.0538 mL	0.2688 mL	0.5376 mL
50 mM	0.0108 mL	0.0538 mL	0.1075 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481