

Coomassie Blue Fast Staining Kit

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at RT

Actual storage temperature shall be subject to the COA.

Biological Description

Description

The Coomassie Blue Fast Staining Kit is used for routine staining and decolorization of protein electrophoresis gels such as SDS-PAGE or non-denaturated PAGE, or for the detection of residual proteins on PAGE gels after Western blotting. Protein dyes based on Coomassie staining are simple and easy to use, allowing for the visualization of proteins separated by electrophoresis with ease.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481