


Tetrahydrojervine

Chemical Properties

CAS No. :	24508-94-9
Formula:	C ₂₇ H ₄₃ NO ₃
Molecular Weight:	429.64
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tetrahydrojervine is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3275 mL	11.6377 mL	23.2753 mL
5 mM	0.4655 mL	2.3275 mL	4.6551 mL
10 mM	0.2328 mL	1.1638 mL	2.3275 mL
50 mM	0.0466 mL	0.2328 mL	0.4655 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Martin HJ, Pachaly P, Zymalkowski F. [Isolation and structural elucidation of quaternary tetrahydroberine alkaloids (author's transl)]. Arch Pharm (Weinheim). 1977 Apr;26(4):314-9. German. PubMed PMID: 860900.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481