

(R)-2-Hexanol

Chemical Properties


CAS No. : 26549-24-6

Formula: C₆H₁₄O

Molecular Weight: 102.18

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(R)-2-Hexanol is a chiral alcohol suitable for biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	9.7867 mL	48.9333 mL	97.8665 mL
5 mM	1.9573 mL	9.7867 mL	19.5733 mL
10 mM	0.9787 mL	4.8933 mL	9.7867 mL
50 mM	0.1957 mL	0.9787 mL	1.9573 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Zhou F, et al. Bacterial volatile ammonia regulates the consumption sequence of d-pinitol and d-glucose in a fungus associated with an invasive bark beetle. ISME J. 2017 Dec;11(12):2809-2820.

Kondo K, et al. Characterization of orally active nonpeptide vasopressin V(2) receptor agonist. Synthesis and biological evaluation of both the (5R)- and (5S)-enantiomers of 2-[1-(2-Chloro-4-pyrrolidin-1-yl-benzoyl)-2,3,4,5-tetrahydro-1H-1-benzazepin-5-yl]-N-isopropylacetamide. J Med Chem. 2002 Aug 15;45(17):3805-8.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481