


5-Methyl-1,2,3,4-tetrahydronaphthalene

Chemical Properties

CAS No. :	2809-64-5
Formula:	C ₁₁ H ₁₄
Molecular Weight:	146.23
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-Methyl-1,2,3,4-tetrahydronaphthalene is a naphthalene derivative suitable for biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.8385 mL	34.1927 mL	68.3854 mL
5 mM	1.3677 mL	6.8385 mL	13.6771 mL
10 mM	0.6839 mL	3.4193 mL	6.8385 mL
50 mM	0.1368 mL	0.6839 mL	1.3677 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481