

1× DPBS Buffer (Sterile)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at 4°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	1× DPBS Buffer (Sterile), also known as Dulbecco's phosphate-buffered saline, is a commonly used buffered saline solution that maintains extracellular fluid osmolarity and pH stability.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481