


2-Ethynyl-6-hydroxy-3-iodopyridine

Chemical Properties

CAS No. :	1824056-39-4
Formula:	C ₇ H ₄ INO
Molecular Weight:	245.02
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Intermediates and Building Blocks - Nucleoside base, Electrophile; Heterocyclic Compounds - Pyridine; Scaffolds and Templates
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0813 mL	20.4065 mL	40.813 mL
5 mM	0.8163 mL	4.0813 mL	8.1626 mL
10 mM	0.4081 mL	2.0406 mL	4.0813 mL
50 mM	0.0816 mL	0.4081 mL	0.8163 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481