


CPI-455 analogue

Chemical Properties

CAS No. :	2309755-73-3
Formula:	C ₁₆ H ₁₄ N ₄ O
Molecular Weight:	278.31
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	CPI-455 analogue (CPI-455 analogue) is an analogue of CPI-455 which is a KDM5 inhibitor.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: insoluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5931 mL	17.9656 mL	35.9312 mL
5 mM	0.7186 mL	3.5931 mL	7.1862 mL
10 mM	0.3593 mL	1.7966 mL	3.5931 mL
50 mM	0.0719 mL	0.3593 mL	0.7186 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Vinogradova M., et al. An inhibitor of KDM5 demethylases reduces survival of drug-tolerant cancer cells. Nat Chem Biol. 2016 Jul;12(7):531-8.
Benjamin R. Leadem. NOVEL HISTONE DEMETHYLASE INHIBITORS SYNERGISTICALLY

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481