


2-Amino-3'-O-methyladenosine

Chemical Properties

CAS No. : 80791-88-4
 Formula: C₁₁H₁₆N₆O₄
 Molecular Weight: 296.28
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Amino-3'-O-methyladenosine is a Nucleoside Derivative - 3'-Modified nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3752 mL	16.8759 mL	33.7519 mL
5 mM	0.675 mL	3.3752 mL	6.7504 mL
10 mM	0.3375 mL	1.6876 mL	3.3752 mL
50 mM	0.0675 mL	0.3375 mL	0.675 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481