

VER-00158411

Chemical Properties


CAS No. : 1174664-88-0

Formula: C₃₁H₃₄N₆O₃

Molecular Weight: 538.64

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	VER-00158411 is a checkpoint kinase 1 and CHK2 inhibitor [IC ₅₀ : 4.4 nM and 4.5 nM, respectively].
Targets(IC ₅₀)	Chk

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8565 mL	9.2826 mL	18.5653 mL
5 mM	0.3713 mL	1.8565 mL	3.7131 mL
10 mM	0.1857 mL	0.9283 mL	1.8565 mL
50 mM	0.0371 mL	0.1857 mL	0.3713 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481