

## EZH2/KMT6 Protein, Human, Recombinant (GST)

### General Information

Protein Construction: Full length of human EZH2 with an N-terminal GST tag.

Species: Human

Expression Host: Baculovirus-Insect Cells

Accession: Q15910

Molecular Weight: 130 kDa (reducing conditions)

### QC Testing

Formulation: Supplied as sterile 50 mM Tris-HCl, pH 7.5, 50-300 mM NaCl, 10 mM glutathione, 0.1 mM EDTA, 0.25-1.0 mM DTT, 0-0.1 mM PMSF, 10-25% glycerol

### Preparation and Storage

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

Proteins are shipped with blue ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481