

Ac-IEVDIDV (TFA)

Chemical Properties

CAS No. :

Formula: C39H62F3N7O17

Molecular Weight: 957.94

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Ac-IEVDIDV TFA is a short peptide sequence.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0439 mL	5.2195 mL	10.4391 mL
5 mM	0.2088 mL	1.0439 mL	2.0878 mL
10 mM	0.1044 mL	0.522 mL	1.0439 mL
50 mM	0.0209 mL	0.1044 mL	0.2088 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481