

Ac-IEVDID (TFA)

Chemical Properties

CAS No. :

Formula: C₃₄H₅₃F₃N₆O₁₆

Molecular Weight: 858.81

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Ac-IEVDID TFA is a short peptide sequence with Ac at the end.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1644 mL	5.822 mL	11.644 mL
5 mM	0.2329 mL	1.1644 mL	2.3288 mL
10 mM	0.1164 mL	0.5822 mL	1.1644 mL
50 mM	0.0233 mL	0.1164 mL	0.2329 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481