


## Nelutroctiv

## Chemical Properties

CAS No. :	2299177-09-4
Formula:	C <sub>24</sub> H <sub>22</sub> F <sub>5</sub> N <sub>3</sub> O <sub>4</sub> S
Molecular Weight:	543.51
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Nelutroctiv is a potent activator of cardiac troponin.
Targets(IC50)	Others, Calcium Channel

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8399 mL	9.1995 mL	18.3989 mL
5 mM	0.368 mL	1.8399 mL	3.6798 mL
10 mM	0.184 mL	0.9199 mL	1.8399 mL
50 mM	0.0368 mL	0.184 mL	0.368 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481