


Compound FL0050

Chemical Properties

CAS No. : 1060816-55-8
 Formula: C₇H₄ClN₃O₂
 Molecular Weight: 197.58
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound FL0050, with CAS No. 1060816-55-8, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Compound FL0050 provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0612 mL	25.3062 mL	50.6124 mL
5 mM	1.0122 mL	5.0612 mL	10.1225 mL
10 mM	0.5061 mL	2.5306 mL	5.0612 mL
50 mM	0.1012 mL	0.5061 mL	1.0122 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481