

Beta-Catenin Protein, Human, Recombinant (aa 2-781, His)

General Information

Synonyms: MRD19;CTNNB;beta-catenin;armadillo;EVR7;NEDSDV;PRO2286;OK/SW-cl.35

Protein Construction: 2-781 aa

Species: Human

Expression Host: E. coli

Accession: P35222

Molecular Weight: 90.9 kDa (Predicted)

AA Sequence:

```
ATQADLMELDMAMEPDRKAAVSHWQQQSYLDSGIHSGATTTAPSLSGKGNPEEEDVDTSQVLYEWEQGFS
QSFTQEQVADIDGQYAMTRAQRVRAAMFPETLDEGMQIPSTQFDAAHPTNVQRLAEPSQMLKHAVVNLIN
YQDDAELATRAIPELTKLLNDEDQVVVNKAAVMVHQLSKKEASRHAIMRSPQMVSIAIVRTMQNTNDVETAR
CTAGTLHNLSSHREGLLAIFKSGGIPALVKMLGSPVDSVLFYAITTLHNLLHQEGAKMAVRLAGGLQKMVAL
LNKTNVKFLAITTDCLQILAYGNQESKLIILASGGPQALVNIMRTYTYEKLLWTTSRVLKVLSSCNKPAIVEAG
GMQALGLHLTDPSQRLVQNCLWTLRNLSDAATKQEGMEGLLGLVQLLGSDDINVTCAAGILSNLTCNNYK
NKMMVCQVGGIEALVRTVLRAGDREDITEPAICALRHLSRHQEAEMAQNAVRLHYGLPVVVKLLHPPSHWP
LIKATVGLIRNLALCPANHAPLREQGAIPRLVQLLVRHQDTQRRTSMGGTQQQFVEGVRMEEIVEGCTGALH
ILARDVHNRIVIRGLNTIPLFVQLLYSPIENIQRVAAGVLCELAQDKEAAEAIEAEGATAPLTELLHSRNEGVATY
AAAVLFRMSEDKPQDYKKRLSVELTSSLFRTEPMAWNETADLGLDIGAQGEPLGYRQDDPSYRSFHSGGYG
QDALGMDPMMEHEMGGHHPGADYPVDGLPDLGHAQDLMDGLPPGDSNQLAWFDTDL
```

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it.

Purity: > 85% as determined by SDS-PAGE

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation:

If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.

If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481