

TROP-2 Protein, Human, Recombinant (His) V2

General Information

Synonyms:	TACD2;EGP1;tumor-associated calcium signal transducer 2;GA733-1;TROP-2;GP50;M1S1;EGP-1;TROP2;GA7331
Protein Construction:	His27-Thr274
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P09758
Molecular Weight:	28.9 kDa (Predicted); 40-60 kDa (Due to glycosylation)

QC Testing

Biological Activity:	Immobilized Human TROP-2, His Tag at 0.5 µg/ml (100 µl/Well). Dose response curve for Anti-TROP-2 Antibody, hFc Tag with the EC50 of 2.2 ng/ml determined by ELISA (QC Test).
Purity:	> 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 µg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Trop-2, also known as epithelial glycoprotein-1 antigen (EGP-1), is a protein that in humans is encoded by the TACSTD2 gene. Mutations of this gene result in gelatinous drop-like corneal dystrophy, an autosomal recessive disorder characterized by severe corneal amyloidosis leading to blindness.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481