


Palmitoyl tripeptide-5 biTFA

Chemical Properties

CAS No. :	623172-56-5
Formula:	C37H67F6N5O9
Molecular Weight:	839.96
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Palmitoyl tripeptide-5 biTFA is a bioactive peptide with anti-aging properties, reported for use as a cosmetic ingredient.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1905 mL	5.9527 mL	11.9053 mL
5 mM	0.2381 mL	1.1905 mL	2.3811 mL
10 mM	0.1191 mL	0.5953 mL	1.1905 mL
50 mM	0.0238 mL	0.1191 mL	0.2381 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481