

Cy3.5 acetate

Chemical Properties

CAS No. :

Formula: C₄₁H₄₆N₂O₁₆S₄

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3.5 acetate is a fluorescent dye used for labeling antibodies.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481