

Sulfo-Cy3(Me)COOH TEA

Chemical Properties

CAS No. :

Formula: C₃₆H₅₁N₃O₈S₂

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sulfo-Cy3(Me)COOH TEA is a derivative of sulfo-Cyanine3. Cyanine3 serves as an orange fluorescent labeling dye for proteins and nucleic acids, with an excitation wavelength (Ex) of 554 nm and an emission wavelength (Em) of 568 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481