

N3-PEG4-YPYDVPDYA-Doxorubicin

Chemical Properties

CAS No. :

Formula: C40H54D2N4O16P2VY3

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	N3-PEG4-YPYDVPDYA-Doxorubicin forms part of an antibody-drug conjugate, consisting of the anthracycline antibiotic Doxorubicin and a degradable ADC linker, N3-PEG4-YPYDVPDYA.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481