

## Pomalidomide-amido-C3-piperazine-N-Boc

### Chemical Properties

CAS No. :

Formula: C<sub>27</sub>H<sub>33</sub>N<sub>5</sub>O<sub>8</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Pomalidomide-amido-C3-piperazine-N-Boc is a synthetic E3 ligase ligand-linker conjugate that includes a cereblon ligand based on Pomalidomide and a linker. It is utilized in the synthesis of PROTAC PD-1/PD-L1 degrader-1.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481