

## Benzyl-PEG8-NHS ester

### Chemical Properties

CAS No. :

Formula: C<sub>28</sub>H<sub>43</sub>N<sub>0</sub>O<sub>12</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Benzyl-PEG8-NHS ester is a PROTAC linker belonging to the PEG class, utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481