

AURKA-IN-1

## Chemical Properties

CAS No. :

Formula: C<sub>25</sub>H<sub>25</sub>FN<sub>8</sub>O<sub>5</sub>S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	AURKA-IN-1 (compound 9) is a cellularly active covalent inhibitor of AURKA and serves as a lysine-reactive activity probe.
Targets(IC50)	Aurora Kinase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481