

N3-PEG4-DYKDDDD-Doxorubicin

Chemical Properties

CAS No. :

Formula: C40H54D5KN4O16Y

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	N3-PEG4-DYKDDDD-Doxorubicin is a component of an antibody-drug conjugate, comprising the anthracycline antibiotic doxorubicin and the degradable ADC linker N3-PEG4-DYKDDDD.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481