

Fmoc-NH-PEG5-NH-Boc

Chemical Properties

CAS No. :

Formula: C₃₂H₄₆N₂O₉

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fmoc-NH-PEG5-NH-Boc is a PROTAC linker that belongs to the PEG class, utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481