

Ms-PEG3-NHS ester

## Chemical Properties

CAS No. :

Formula: C<sub>14</sub>H<sub>23</sub>N<sub>1</sub>O<sub>10</sub>S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Ms-PEG3-NHS ester is a PROTAC linker belonging to the PEG class and is suitable for synthesizing PROTAC molecules.
Targets(IC50)	PROTAC Linker

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481