

Benzyl-PEG24-azide

Chemical Properties

CAS No. :

Formula: C55H103N3O24

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Benzyl-PEG24-azide is a PROTAC linker within the PEG category, used for synthesizing PROTAC molecules. As a click chemistry reagent, it features an azide group that participates in copper-catalyzed azide-alkyne cycloaddition reactions (CuAAC) with alkyne-containing molecules. Additionally, it engages in strain-promoted azide-alkyne cycloaddition reactions (SPAAC) with molecules having DBCO or BCN groups.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481