

Boc-NH-PEG1-C5-OH

Chemical Properties

CAS No. :

Formula: C₁₂H₂₅N₁O₄

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Boc-NH-PEG1-C5-OH is a PROTAC linker belonging to the PEG category, used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481