

N-Ethyl-3,3,3-trifluoro-N-methylpropanamide-PEG2-Br

Chemical Properties

CAS No. :

Formula: C₁₀H₁₇BrF₃NO₃

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	N-Ethyl-333-trifluoro-N-methylpropanamide-PEG2-Br is a PEG-based PROTAC linker, utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481