

## Thalidomide-NH-C5-NH2

## Chemical Properties

CAS No. : 2093386-39-9


Formula: C18H22N4O4

Molecular Weight: 358.39

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Thalidomide-NH-C5-NH2 is a synthetic E3 ligase ligand-linker conjugate, consisting of a Thalidomide-based cereblon ligand and a linker, used in the synthesis of PROTACs.
Targets(IC50)	Apoptosis, Autophagy, E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7903 mL	13.9513 mL	27.9026 mL
5 mM	0.5581 mL	2.7903 mL	5.5805 mL
10 mM	0.279 mL	1.3951 mL	2.7903 mL
50 mM	0.0558 mL	0.279 mL	0.5581 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481