

## Propargyl-PEG9-THP

### Chemical Properties

CAS No. :

Formula: C<sub>26</sub>H<sub>48</sub>O<sub>11</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Propargyl-PEG9-THP is a PROTAC linker in the PEG class, useful for the synthesis of PROTAC molecules. This compound is a click chemistry reagent containing an alkyne group, enabling it to undergo copper-catalyzed azide-alkyne cycloaddition reactions (CuAAC) with azide-containing molecules.
Targets(IC50)	PROTAC Linker

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481