

## Biotin-PEG3-methyl ethanethioate

### Chemical Properties

CAS No. :

Formula: C<sub>20</sub>H<sub>35</sub>N<sub>3</sub>O<sub>6</sub>S<sub>2</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Biotin-PEG3-methyl ethanethioate is a PROTAC linker that belongs to the PEG category and is used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481