

Thiol-PEG-COOH (MW 5000)

Chemical Properties

CAS No. :

Formula: C₁₁H₂₁N₀S₅

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Thiol-PEG-COOH (MW 5000) is a PROTAC linker within the PEG category. Boc-NH-C6-amido-C4-acid is employed in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481