

Fmoc-N-amido-PEG36-Boc

Chemical Properties

CAS No. :

Formula: C₉₄H₁₆₉N₄O₄₀

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fmoc-N-amido-PEG36-Boc is a PROTAC linker that belongs to the PEG class and is used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481