

## Benzyl-PEG11-Boc

### Chemical Properties

CAS No. :

Formula: C<sub>35</sub>H<sub>62</sub>O<sub>14</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Benzyl-PEG11-Boc is a PROTAC linker within the PEG category. It is utilized for the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481