

Boc-NH-PEG11-C2-acid

Chemical Properties

CAS No. :

Formula: C₃₀H₅₉N₁₅O₁₅

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Boc-NH-PEG11-C2-acid is a PEG-based PROTAC linker, utilized in the synthesis of PROTAC molecules. |
| Targets(IC50) | PROTAC Linker |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481