

## m-PEG7-Silane

## Chemical Properties

CAS No. :

Formula: C<sub>25</sub>H<sub>53</sub>NO<sub>11</sub>Si

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	m-PEG7-Silane is a PROTAC linker belonging to the PEG category and is used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481