

Fmoc-NH-PEG14-acid

## Chemical Properties

CAS No. :

Formula: C46H73NO18

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | Fmoc-NH-PEG14-acid is a type of PROTAC linker within the PEG class, and it is used in the synthesis of PROTAC molecules. |
| Targets(IC50) | PROTAC Linker  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481